

MASA BAĞLANTILI EKARTÖRLER
TABLE CONNECTION RETRACTORS

Jinekoloji / GYNAECOLOGY - **JN**





L 630.22
Bağlantı klemp/Connecting clamp



L 630.24
Çentik bağlantı çubuğu
Extension for notch rod
10 cm



ROCHARD
L 630.29
Adaptör/Adapter



L 630.23
Çentik Çubuğu/Notch rod
15 cm



L 630.25
Çentik son bağlantı parçası
Notch rod end part
6 cm



L 630.27
Düz bağlantı çubuğu
Extension rod, smooth
20 cm



K 742.00

Ameliyat masalarının her iki yan demirlerine bağlama için montaj çerçevesi. Daha iyi sağlamlık sağlanmış olmakta olup sterilizasyon için sökülebilir.
Mounting frame for fastening to both sides of the lateral bars of the operating tables, thus more stability is given, can be dismantled for sterilization.



DEAVER
L 630.30
Ekartör/Retractor
85 x 25 mm



MIKULICZ
L 630.31 120 x 50 mm
L 630.32 155 x 50 mm
Ekartör/Retractor



KELLY
L 630.34
Ekartör/Retractor
140 x 40 mm



HARRINGTON
L 630.40
Ekartör/Retractor
130 x 62 mm



L 630.37
Karaciğer ekartörü/Liver Retractor
140 x 85 mm



L 630.38
İdrar kesesi spatulası
Bladder spatula
30° açılı/30° angled
110 x 85 mm

