


**MİKRO AÇILI DAMAR MAKASLARI**  
MICRO ANGLED VASCULAR SCISSORS

Kalp Damar ve Göğüs Cerrahi/ Cardiovascular and Thoracic Surgery - KR



Yuvarlak Sap  
Round Handle

Ince Bleyd  
Fine Blade



Hassas Bleyd, tek tarafı küre uçlu bleyd  
Delicate Blade, with one probe pointed blade



Hassas Bleyd  
Delicate Blade



Ultra Hassas Bleyd  
Ultra Delicate Blade



VI

C 169.00 - C 169.56  
16.5 cm