

EL MATKAPLARI
HAND DRILLS

Ortopedi / Orthopedics - **OR**





V 606.02
Anahtar, bilye kilitli
Key, with ball lock

BUNNEL
V 615.00


Küçük el matkabi, KIRSCHNER tel için delikli küçük el matkabi,
3 çene ağızlı, 4 mm çaplı başlıklar için uygundur.
Small hand drill, canulated for KIRSCHNER wires, with 3-jaw chuck,
to fit any size of shanks up to 4.0 mm diameter.

V 606.02 anahtarı ile
With key V 606.02
16 cm



Ayarlı sap
Adjustable handle

V 612.00
El matkabi/Hand drill only
28 cm



V 606.01
Anahtar/Key

El matkabi paslanmaz, 3 çene ağızlı, 6.35 mm çaplı başlıklar için uygundur. Orjinal 'Jacops Chuck' aynaları ile ve V 606.01 anahtarı ile, anahtarı veya anahtarsız kullanılabilir.
Hand drill, stainless, original 'Jacops Chuck' and key V 606.01
Can be used with or without the key.