

**ASPIRATÖR UÇLARI**  
SUCTION TUBES

Plastik Cerrahi / Plastic Surgery - PC





**BARON**  
**J 234.00** Ø 1.5 mm  
7.5 cm

**BARON**  
**J 234.01** Ø 2.0 mm  
7.5 cm

**BARON**  
**J 234.02** Ø 1.0 mm  
7.5 cm

**SCHUKNECHT**  
**J 370.00** Ø 1 mm  
**J 370.01** Ø 1.5 mm  
**J 370.02** Ø 2 mm  
**J 370.03** Ø 3 mm  
13 cm

Ø 5.0 mm



**J 230.00**  
LUER  
RECORD

**J 231.00**  
RECORD  
LUER-LOCK

Olive  
Ø 6-9 mm

1/2

**J 352.00**  
Uc genişliği 6.5 mm  
Tip wide 6.5 mm  
LUER bağlantılı  
With LUER hub  
18 cm



Olive  
Ø 6-9 mm

1/2

**FERGUSSON**  
**J 220.02** Ø 4 mm  
İzoleli/Insulated  
LUER bağlantılı  
With LUER hub  
18 cm



Olive  
Ø 6-9 mm

1/2

**FRAZIER**  
**J 222.00** Ø 3 mm  
**J 222.01** Ø 5 mm  
LUER bağlantılı  
With LUER hub  
18 cm