


ASPIRATÖR UÇLARI
SUCTION TUBES

Veteriner/ Veterinary - VT





	FERGUSSON				
Çalışma boyu Working length	11 cm	13 cm	16 cm	20 cm	
Tam boy Total Length	18 cm	21.5 cm	24.5 cm	28 cm	Ø
	J 215.01	J 215.21	J 215.51	J 215.81	1.5 mm
	J 215.02	J 215.22	J 215.52	J 215.82	2.0 mm
	J 215.03	J 215.23	J 215.53	J 215.83	2.5 mm
	J 215.04	J 215.24	J 215.54	J 215.84	3.0 mm
	J 215.06	J 215.26	J 215.56	J 215.86	4.0 mm
	J 215.08	J 215.28	J 215.58	J 215.88	5.0 mm

Parmak kumandalı, temizleme telli
 parmak basma yüzeyi yatay LUER bağlantı
 With finger cut-of and styled
 Thumb piece 90° to curve, With LUER hub

 Olive
 Ø 6-9 mm


FERGUSSON
J 220.02 Ø 4 mm
 İzoleli, parmak kumandalı,
 temizleme telli
 parmak basma yüzeyi yatay
 LUER bağlantı
 Isolated, with finger cut-of and stylet,
 thumb piece 90° to curve, with LUER
 hub
 18 cm

 Olive
 Ø 6-9 mm


FRAZIER
J 222.00 Ø 3 mm
J 222.01 Ø 5 mm
 Parmak kumandalı, temizleme
 telli parmak basma yüzeyi
 yatay LUER bağlantı
 With finger cut-of and stylet
 Thumb piece 90° to curve,
 with LUER hub
 18 cm

FRAZIER
J 221.00 Ø 2.3 mm
J 221.01 Ø 2.7 mm
J 221.02 Ø 3.4 mm
J 221.03 Ø 4.0 mm
 Parmak kumandalı, temizleme
 telli parmak basma yüzeyi yatay
 LUER bağlantı
 With finger cut-of and stylet
 Thumb piece 90° to curve,
 with LUER hub
 11 cm

